


# Introduction to Fourier Transform Infrared Spectroscopy (FTIR)

FTIR measures infrared absorption and emission to identify chemical bonds.

This non-destructive technique is vital in materials science, chemistry, and pharmaceuticals.





# Fundamentals and Principles of FTIR

## Infrared absorption

Molecules absorb IR radiation at frequencies matching vibration modes.

## Fourier Transform

Converts raw interferogram data into a readable spectrum.

## Vibrational frequencies

Specific to molecular bonds, unique fingerprints for identification.

# Uses and Applications of FTIR


## Material Analysis

- **Polymer characterization**
- **Chemical compound ID**
- **Coating analysis**

## Pharmaceuticals & Food

- **Drug purity checks**
- **Adulteration detection in foods**
- **Environmental pollutant analysis**





# Sample Requirements and Applicable Materials

## Sample Forms

Solids, liquids, and gases can be analyzed.

## Preparation Methods

KBr pellets, thin films, ATR technique.

## Material Types

Polymers, organics, inorganics, coatings.

## Limitations

Strong water absorption and particle size effects.

# FTIR Data Acquisition and Results

## Data Form

Interferogram converted to Absorbance vs. Wavenumber spectrum.

## Peak Analysis

Peaks correspond to molecular vibrations showing chemical groups.

## Preprocessing

Baseline correction and normalization improve accuracy.

# Data Interpretation and Property Calculation

## Identify Functional Groups

Peak positions reveal groups like C=O or N-H.

## Quantify Concentrations

Peak intensity proportional to compound amount.

## Advanced Analysis

Deconvolution resolves overlapping peaks for clarity.

## Calculations

Crystallinity index from peak ratios in polymers.



# Research Example: Polymer Blend

## 1 Study Insight


Peak shifts indicate polymer miscibility.

## 2 Calculations

Blending degree estimated using peak area ratios.

## 3 Conclusion

FTIR confirms molecular interactions in blends.







# Summary and Conclusion



## Powerful Tool

Essential for diverse material characterization.



## Expansive Applications

Used in polymers, pharmaceuticals, environment, and food.



## Data Interpretation

Requires understanding of vibrational spectra and calculations.



## Advancements

Ongoing improvements enhance accuracy and sensitivity.